


MAX VPC™ Screw System

Product Brochure




A big step forward for small bone fractures.



The MAX VPC (Variable Pitch Compression) Screw System is comprised of headless, fully threaded, cannulated screws designed with a constant pitch thread form near the tip and a variable pitch thread form near the head of the screw giving the user tactile feel during installation, all while providing compression to the bone fragments.

The tip of the screw features four forward cutting flutes and two side cutting flutes designed to provide self-tapping capabilities.



Cutting Flutes

Side Cutting

Forward Cutting




The smaller diameter head size requires reduced bone removal during insertion.

The fully threaded design provides bone purchase throughout the entire length of the screw which results in structural support during bone healing.

The screws are driven with a standard metric Hex drive (1.5 mm, 2.0 mm, 2.5 mm) that incorporates a tapered tip to secure the screw to the driver.

Cobalt chrome K-wires can provide increased stiffness and reduces the likelihood of kinking during use.*

The shorter lengths (95 mm, 105 mm, 127 mm) make instrumentation more user friendly and facilitate reduction by securely holding bone fragments.



The instrumentation was developed around the simplicity of the surgery. The color coded system allows for easy identification of specialized instrumentation to the corresponding sized screws.

* Yield strength for Stainless Steel is 110 Ksi min versus 208 Ksi for 35NLT Alloy CoCr.

INDICATIONS

The Zimmer Biomet Variable Pitch Compression Screw System is indicated for alignment and stabilization of small bone fractures.

Specifically:

- Fixation of small bones, such as those in the foot, ankle, wrist, elbow and hand for treatment of fractures, non-unions, or mal-unions
- Ligament reconstruction
- Osteochondritis dissecans
- Arthrodesis of the foot, ankle, wrist, elbow and hand
- Small bone osteotomies, including first metatarsal head osteotomy, metatarsal osteotomies, phalangeal osteotomies, and carpal/metacarpal osteotomies.

All content herein is protected by copyright, trademarks, and other intellectual property rights owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

This material is intended for health care professionals and the Zimmer Biomet sales force. Distribution to any other recipient is prohibited. For complete product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see the package insert and www.zimmerbiomet.com

Check for country product clearances and reference product specific instructions for use. Not for distribution in France.

©2016 Zimmer Biomet



Authorized Representative

Zimmer Biomet UK Ltd.
Waterton Industrial Estate
Bridgend, South Wales
CF31 3XA
UK



ZIMMER BIOMET

Your progress. Our promise.™



Legal Manufacturer

Zimmer Biomet
P.O. Box 587
56 East Bell Drive
Warsaw, Indiana 46581-0587
USA